HT series Compact Scale Operating Instructions

Cautions

- Avoid installing the scale in direct sunlight.
- Avoid over / shock load that could cause damage to the scale.
- Use the scale under conditions of -10°C ~ 40°C, less than 85% R.H.
- Do not mix the type of batteries. Replace the batteries all together.
- If the scale is not to be used for a long period of time, remove all batteries.
- Use only the specified AC adapter for the scale (AC adapter is optional).

Specifications

| | |
|--------------|--|
| HT-300 | 310 x 0.1 g / 10.9 x 0.01 oz |
| HT-3000 | 3100 x 1 g / 109 x 0.1 oz / 6.8 lb x 0.1 oz |
| HT-500 | 510 x 0.1 g / 17.9 x 0.01 oz / 1.1 lb x 0.01 oz |
| HT-5000 | 5100 x 1 g / 179 x 0.1 oz / 11 lb x 0.1 oz |
| Power | 4 x "AA" size batteries or AC adapter (option) |
| Battery life | Approximately 450 hours with alkaline cells at 20°C / 68°F |
| Pan size | 132 (W) x 130 (D) mm / 5.2 (W) x 5.1 (D) in. |
| Dimensions | 136 (W) x 195 (D) x 47 (H)mm / 5.4 (W) x 7.7 (D) x 1.9 (H) in. |

Keys and Symbols

| | |
|--|---|
| ON/OFF key | Turns the scale power on/off. |
| RE-ZERO key | Sets the display to zero. |
| UNITS key | Changes weighing units, g, oz or lb-oz. |
| Symbol | Indicates when the reading is stable. |
| Symbol ZERO ◀ | Indicates when the scale zero is correct. |
| Symbol NET ◀ | Indicates when the net weight is displayed. |
| Symbol | Battery remaining level, shown after power-on only. |

Installing the batteries / Connecting the AC adapter


Installing / Replacing the batteries

Remove the battery cover. Insert 4 new batteries (AA size) into the battery compartment, taking extreme care of the polarities. Replace the battery cover.

The batteries provided with the scale are for testing the performance of the scale and may have a limited life.

AC adapter (optional)

Verify that the AC adapter is correct. Open the rubber cap and plug the AC adapter into the AC adapter jack.

Weighing procedures

1. Turning the scale on / off

Press the **ON/OFF** key to turn the scale on. The scale displays all segments for a few seconds and then displays zero. Press the **ON/OFF** key again to turn the scale off.

2. Selecting the weighing unit

Press the **UNITS** key to select the weighing unit. The scale will power up with the weighing unit that was in use when turned off last time.

3. Weighing

Verify that the reading is zero. If not, press the **RE-ZERO** key to zero the display. If you use a container for weighing, place an empty container on the weighing pan and press the **RE-ZERO** key to zero the display. Place the objects to be weighed on the pan or in the container. When the reading becomes stable, the STABLE indicator is displayed.

Note

Automatic power-off function

If the scale is left ON and the STABLE indicator is displayed, the automatic power-off function turns power off after approximately 5 minutes. To disable or enable this function, press the **ON/OFF** key while pressing the **RE-ZERO** key. Then **POFF** (disabled) or **PoFF** (enabled) will be displayed. Select the setting by pressing the **RE-ZERO** key and press the **UNITS** key. The scale displays **End** and returns to the weighing mode.

Error message

| | |
|---|---|
| F | Overload. Remove the objects from the pan. |
| ---- | The scale is out of zero range or unstable when powered on. |
| Lb | Low battery. Replace the batteries or use the AC adapter. |

Calibration and other functions

Please visit the A&D website (<http://www.aandd.jp/>).

A&D

A&D Company, Limited
3-23-14 Higashi-Ikebukuro,
Toshima-ku, Tokyo 170-0013 JAPAN
<http://www.aandd.co.jp/>

A&D AUSTRALASIA PTY LTD
32 Dew Street, Thebarton, South
Australia 5031 AUSTRALIA
<http://www.andweighing.com.au/>

1WMFPA4001904C

A&D ENGINEERING, INC.
1756 Automation Parkway, San Jose,
California 95131 U.S.A.
<http://www.andonline.com/>

A&D KOREA Limited
817, Manhattan Bldg., 33,
Gukjegeumyung-ro 6-gil,
Yeongdeungpo-gu, Seoul, KOREA
<http://www.andk.co.kr/>

A&D INSTRUMENTS LIMITED
Unit 24/26 Blacklands Way, Abingdon
Business Park, Abingdon, Oxfordshire
OX14 1DY United Kingdom
<http://www.aandd-eu.net/>

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A&D ENGINEERING, INC.
1756 Automation Parkway, San Jose, California 95131 U.S.A.
http://www.andonline.com/

A&D AUSTRALASIA PTY LTD
32 Dew Street, Thebarton, South Australia 5031 AUSTRALIA
http://www.andweighing.com.au/

HT series Compact Scale

Calibration and Function Setting Procedures

Calibration

When is Calibration Required?

Calibration may be required when the HT series compact scale was initially installed or has been moved to another location.

Calibration Using a Calibration Weight


1. Entering the calibration mode

Remove the calibration switch cover located on the bottom of the scale.

Press the [ON:OFF] switch to turn the scale on.

Press the calibration switch while the scale is in the weighing mode.

`[CAL]` will be displayed.

2. Zero calibration

Press the [RE-ZERO] key. `[CALD]` will be displayed.

Wait for the STABLE Indicator (O) to be displayed. Press the [RE-ZERO] key to perform zero calibration.

`[CALF]` will be displayed after a few seconds. To calibrate span, go to step 3.

To return to the weighing mode without performing span calibration, press the [UNITS] key.

3. Span calibration

When `[CALF]` is displayed, place the calibration weight on the center of the weighing pan. Wait for the STABLE indicator to be displayed. Press the [RE-ZERO] key to perform span calibration.

`[End]` will be displayed and the scale will automatically return to the weighing mode.

Note: For details about the calibration weight, see "Specifications".

Calibration by Gravity Compensation

If the acceleration of gravity at your location is not 9.798 m/s^2 and you do not have calibration weights, the scale can be calibrated by compensating for the acceleration of gravity. (See "The Value of Gravity at Various Locations".)

1. Setting a new acceleration value

When `[CAL]` is displayed, press the [UNITS] key. `[9798]` will be displayed.

To change the value, press the [RE-ZERO] key to increment the blinking digit and press the [UNITS] key to move the blinking digit.

2. Storing the value into the memory

While pressing the [UNITS] key, press and hold the [RE-ZERO] key and release the [UNITS] key.


After `[End]`, `[CAL]` will be displayed. Turn the scale off to finish the procedure.

Specifications

| MODEL | HT-300 | HT-3000 | HT-500 | HT-5000 |
|---------------------|---------------|--------------|---------------|--------------|
| Capacity/Resolution | 310 g x 0.1 g | 3100 g x 1 g | 510 g x 0.1 g | 5100 g x 1 g |
| Maximum tare | 310 g | 3100 g | 510 g | 5100 g |
| Calibration weight | 300g ± 0.01g | 3000g ± 0.1g | 500g ± 0.01g | 5000g ± 0.1g |


The Value of Gravity at Various Locations

| | | | | | |
|--------------|------------------------|--------------------|------------------------|----------------|------------------------|
| Amsterdam | 9.813 m/s ² | Havana | 9.788 m/s ² | Rio de Janeiro | 9.788 m/s ² |
| Athens | 9.807 m/s ² | Helsinki | 9.819 m/s ² | Rome | 9.803 m/s ² |
| Auckland NZ | 9.799 m/s ² | Kuwait | 9.793 m/s ² | San Francisco | 9.800 m/s ² |
| Bangkok | 9.783 m/s ² | Lisbon | 9.801 m/s ² | Singapore | 9.781 m/s ² |
| Birmingham | 9.813 m/s ² | London (Greenwich) | 9.812 m/s ² | Stockholm | 9.818 m/s ² |
| Brussels | 9.811 m/s ² | Los Angeles | 9.796 m/s ² | Sydney | 9.797 m/s ² |
| Buenos Aires | 9.797 m/s ² | Madrid | 9.800 m/s ² | Taichung | 9.789 m/s ² |
| Calcutta | 9.788 m/s ² | Manila | 9.784 m/s ² | Taiwan | 9.788 m/s ² |
| Cape Town | 9.796 m/s ² | Melbourne | 9.800 m/s ² | Taipei | 9.790 m/s ² |
| Chicago | 9.803 m/s ² | Mexico City | 9.779 m/s ² | Tokyo | 9.798 m/s ² |
| Copenhagen | 9.815 m/s ² | Milan | 9.806 m/s ² | Vancouver, BC | 9.809 m/s ² |
| Cyprus | 9.797 m/s ² | New York | 9.802 m/s ² | Washington DC | 9.801 m/s ² |
| Djakarta | 9.781 m/s ² | Oslo | 9.819 m/s ² | Wellington NZ | 9.803 m/s ² |
| Frankfurt | 9.810 m/s ² | Ottawa | 9.806 m/s ² | Zurich | 9.807 m/s ² |
| Glasgow | 9.816 m/s ² | Paris | 9.809 m/s ² | | |

Units

All the available units and the order they appear in the display are as follows:

Select either one of Taiwan tael or Hong Kong tael.

The units to be used can be selected in the function setting mode. The displaying order of the selected units is the same as above, while skipping the units that are not selected.

For detailed information about selecting units, see “Function Setting”.

Function Setting

Selecting Units

1. While pressing the [UNITS] key, press the [ON:OFF] key to turn the scale on. **Unit** will be displayed.
2. Press the [RE-ZERO] key once. **Unit g** will be displayed.
3. Each time the [RE-ZERO] key is pressed, the display switches between **Unit g** and **Unit g**.

Note: The STABLE indicator in the upper left corner of the display indicates that the unit is selected and can be used for weighing.

4. Press the [UNITS] key to select a unit. Units are displayed in turn, as shown below:
tola **t** → ounce (avoir) **oz** → troy ounce **ozt** → pound ounce **lb oz** → Taiwan tael **tl**
→Hong Kong tael **tl**.
5. Select a unit to be used and press the [RE-ZERO] key to display the STABLE indicator.
6. Press the [UNITS] key. **End** will be displayed blinking.
Press the [RE-ZERO] key. The new setting will be stored and the display will stop blinking. After a second, the scale will automatically return to the weighing mode.

Changing the Decimal Point Type / Changing the Filter

1. While pressing the [UNITS] key, press the [ON:OFF] key to turn the scale on. **Unit** will be displayed.
2. Press the [UNITS] key. **Pnt.** will be displayed. When the decimal point type is not to be changed, go to step 4.
3. Each time the [RE-ZERO] key is pressed, the display switches between **Pnt .** and **Pnt ,**.
Select the decimal point type to be used.
4. Press the [UNITS] key. **Filter** will be displayed. When the filter is not to be changed, go to step 6.
5. Each time the [RE-ZERO] key is pressed, the display switches between **Filter** and **Filter**.
Select the filter to be used.

Note: **Filter** has been set at the factory before shipping. **Filter** is a filter with different frequency characteristics. Select one appropriate to the operating environment.

6. Press the [UNITS] key. **End** will be displayed blinking.
Press the [RE-ZERO] key. The new setting will be stored and the display will stop blinking. After a second, the scale will automatically return to the weighing mode.