

# TM902

## Animal Scale



## Top Quality, Low Profile Performer from Totalcomp

- *60 x 0.02 lb or 30 x 0.01kg*
- *0.75" LCD display*
- *Keys for Tare, Zero, Hold On/Off*
- *Hold feature for moving animals*
- *Battery or AC operation*
- *Stainless steel platter*
- *Changeable from lb to kg*



**Totalcomp** Scales & Components

13-01 Pollitt Drive Ext, Fair Lawn NJ 07410  
Phone 201-797-2718 Fax 201-797-2287

email: [sales@totalcomp.com](mailto:sales@totalcomp.com)

# Specifications

# TM902

TM902 Animal scale is light weight, and easy to use Veterinary scale. 0.75" LCD display with blue back light. Keys for Tare, Zero, Hold, On/Off. Units are changeable between kg and lb by pressing the "0" button and holding for 2 seconds. Hold feature for animal weighing. While the animal is on the weighing pan, the indicator will hold the weight even if the animal is moving on the platform. Stainless steel removable platter. Adjustable feet with front bubble level. Coil cord between indicator and base is 10' in length. In-use cover included. Powered by internal rechargeable battery. AC adapter for charging included.

Part Number	Capacity	Platform	Weight	Stock
TM902	30 kg x 0.01 kg 60 lb x 0.02 lb	12" x 10"	11 lb	3960

#### **General**

**Weighing resolution**  
1/3,000 divisions

**Accuracy**  
0.01% weighing

**Keyboard**  
Audible key switches

#### **Electrical**

**Power source**  
AC adapter  
100-240VAC  
12VDC, 500ma  
Internal rechargeable  
battery pack (6)AA

**Coil cord**  
10 ft

#### **Temperature**

**Compensated Range**  
-10 deg C to +40 deg C

**Operating Range**  
-16 deg C to +65 deg C

**Weight**  
Net 8.91 lb  
Shipping 10.43 lb

#### **Mechanical**

**Stainless steel platter**  
12" x 10"

**Dimensions**  
Indicator  
8-7/8"Wx2-1/4"Dx5"H  
Base  
12'Lx10"Wx2-3/4"H"

#### **Totalcomp Limited Warranty**

Totalcomp warrants this product to be free of defects in materials and/or workmanship and suitable for the purpose intended as outlined in this sheet. This warranty is effective and shall cover the purchaser for one year from the date of shipment from our plant. If this product is found to be defective by our inspection in accordance with the above listed criteria, we will replace or repair it at our expense. For warranty service, please obtain a return authorization number from us and return the item, shipping prepaid, with a written description of the problem. We will respond promptly with the results of our evaluation.



**Totalcomp Scales & Components**

13-01 Pollitt Drive Ext, Fair Lawn NJ 07410  
Phone 201-797-2718 Fax 201-797-2287

email: [sales@totalcomp.com](mailto:sales@totalcomp.com)