

TJA

Analytical Balance



TJA
Analytical
Balance

Top Quality Performer from Totalcomp

- *Super bright LCD display*
- *Auto calibration*
- *RS232 interface*
- *Stainless steel platter*
- *Breeze break*
- *Units mg, g, oz, ct*

 **Totalcomp** Scales & Components

13-01 Pollitt Drive Ext, Fair Lawn NJ 07410
Phone 201-797-2718 Fax 201-797-2287

email: sales@totalcomp.com

Specifications

TJA

TJA Totalcomp Analytical Balance is ideal for laboratory weighing. Easy to see in low light conditions 9 digit LCD display with backlight. Stainless steel weighing pan 3.15" dia. Sensitivity and speed are separately adjustable in 4 steps from low to high. Countiuous output RS232 via 9 pin connector for remote display or computer input. Below balance weighing for specific gravity determination.

Weighing functions for: Standard weighing and counting.

Keys for: OFF, PRINT, O/T, CAL, ON.

Units are g, mg, ct, oz.

Glass breeze break with top and side doors allow for quick placement of item to be weighed. CAL button and 200g calibration weight allow for easy recalibration of balance. Rear bubble level and two adjustable feet allow for quick leveling. 13.4"L x 8.5"W x 13.8"H. 120 VAC.

Part Number	Capacity	Platform	Weight	Stock
TJA-203	200 g x 1 mg	80mm dia	20 lb	5966

Electrical

Power source

120V AC adapter
12VDC, 500ma

Temperature

Compensated Range

-10 deg C to +40 deg C

Operating Range

-16 deg C to +65 deg C

Dimensions

13.4"L x 8.5"W x 13.8"H

Totalcomp Limited Warranty

Totalcomp warrants this product to be free of defects in materials and/or workmanship and suitable for the purpose intended as outlined in this sheet. This warranty is effective and shall cover the purchaser for one year from the date of shipment from our plant. If this product is found to be defective by our inspection in accordance with the above listed criteria, we will replace or repair it at our expense. For warranty service, please obtain a return authorization number from us and return the item, shipping prepaid, with a written description of the problem. We will respond promptly with the results of our evaluation.



Totalcomp Scales & Components

13-01 Pollitt Drive Ext, Fair Lawn NJ 07410

Phone 201-797-2718 Fax 201-797-2287 email: sales@totalcomp.com